Dutchblue.com Nederland

Markt 30a
2611 GV Delft
Nederland

T 00 800 2611 2611 (gratis) of 015 3030 025
F 084 883 66 06
E info@dutchblue.com
W www.dutchblue.com

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Square de Meeûs 37
1000 Brussel
België

T 00 800 2611 2611 (gratis) of (0032) 026 699 663
F (0032) 02 403 05 55
E info@dutchblue.com
W www.dutchblue.com

Dutchblue.com Deutschland

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Dutchblue.com United Kingdom

Park House
111 Uxbridge Road
Ealing, London, W5 5TL

T 00 800 2611 2611 (free) or 02 033 977 242
F 02 033 977 524
E info@dutchblue.com
W www.dutchblue.com

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User Manual – Hip protector pants with protectors

REF 4000 serie
IFU 4000 – 13M2

dutchblue.com

HIPPIES hip protector pants with protectors

HIPPIES are supplied as a set of four hip protection pants in one size and two protectors or as a set of one hip protection pant with two stitched protectors. The hip protector pants also come as a set of four without protections. The HIPPIES hip protection pants with protectors can be worn comfortably during the day and at night, and can be easily combined with incontinence products.

Composition

Pants	95% cotton 5% lycra
Protectors	D3O® material

Instructions for use

Inserting the protector in the pants

- Both sides of the hip protector pants have been fitted with inside pockets with flaps.
- Before putting on the hip protector pants insert both protectors into the pockets, with the narrow part pointing downwards and the flat side facing inwards.
- The pockets are closed with a flap to ensure that the protectors remain securely in place.
- Once both protectors have been inserted into the pants, the hip protector pants with protectors can be worn in the customary fashion.
- The hip protector pants with sewn-in protectors can be worn immediately in the customary fashion.

Instructions for cleaning

Hip protector pants

The hip protector pants are machine washable at 60°C and can be tumble dried. Always take out the removable protectors before washing the hip protector pants.

Hip protector pant with stitched protectors

The hip protector pants with stitched protectors can be machine washed as a set at 60°C. Do not tumble dry.

Protectors

- Removable protectors can also be cleaned with a soft cloth and mild detergent or clinical cleaning agents.
- Protectors cannot be tumble dried.
- Never use disinfection agents, phenol or solvents.
- Do not use bleaching agents as they can damage the HIPPIES over time.

PLEASE NOTE:

- Regularly check the protectors for damage.
- HIPPIES cannot prevent every hip fracture, but contribute substantially to minimising the risk of contracting a hip fracture.
- Avoid contact with sharp objects.

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