Dutchblue.com Nederland

Markt 30a
2611 GV Delft
Nederland

T 00 800 2611 2611 (gratis) of 015 3030 025
F 084 883 66 06
E info@dutchblue.com
W www.dutchblue.com

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Square de Meeûs 37
1000 Brussel
België

T 00 800 2611 2611 (gratis) of (0032) 026 699 663
F (0032) 02 403 05 55
E info@dutchblue.com
W www.dutchblue.com

Dutchblue.com Deutschland

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Dutchblue.com United Kingdom

Park House
111 Uxbridge Road
Ealing, London, W5 5TL

T 00 800 2611 2611 (free) or 02 033 977 242
F 02 033 977 524
E info@dutchblue.com
W www.dutchblue.com

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User Manual – Anti-decubitus top mattress

REF 3010
IFU 3010 – 13M2

SLEEPIES anti-decubitus top mattress

Sleep peacefully and comfortably on our SLEEPIES anti-decubitus top mattress. The combination of air and foam ensures the perfect distribution of pressure, reducing the possibility of contracting wounds caused by pressure complaints while reducing pain as well. The Dutchblue anti-decubitus top mattress can be used as a preventative measure or therapeutically, depending on the weight and physical condition of the resident. Ideal for the prevention and cure of decubitus complaints.

Instructions for use

Inflating the mattress

- Remove the packaging. Do not use scissors, knives or other sharp implements when opening the packaging.
- Remove the top mattress from the carry bag.
- Remove the bed linen.
- Place the top mattress at the foot of the mattress.
- Open both valves by turning them anti-clockwise three times.
- The top mattress will inflate and unroll automatically until it is positioned flat on the bed mattress.
- The first time this is done the mattress will take about ten minutes to fully inflate. After this, the mattress will be inflated in only 3 to 5 minutes.
- Once the top mattress has been fully inflated, turn the valves clockwise to shut.
- You can then make the bed as usual.

Deflating the mattress

- To deflate the mattress, open the valves completely.
- Roll the mattress up as tightly as possible, working in the direction of the valves.
- You will hear the air coming out of the mattress.
- Once the mattress has been completely rolled up, close the valves securely.
- Store the mattress in the carry bag that was included.

Instructions for cleaning

The mattress can be cleaned with a soft cloth and mild detergent or clinical cleaning products. Never use disinfection agents, phenol or solvents.

PLEASE NOTE:

- The typical odour of visco-elastic foam will disappear in time.
- Do not allow to come into contact with open fire.
- Avoid contact with sharp objects.

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