Dutchblue.com Nederland

Markt 30a
2611 GV Delft
Nederland

T 00 800 2611 2611 (gratis) of 015 3030 025
F 084 883 66 06
E info@dutchblue.com
W www.dutchblue.com

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Square de Meeûs 37
1000 Brussel
België

T 00 800 2611 2611 (gratis) of (0032) 026 699 663
F (0032) 02 403 05 55
E info@dutchblue.com
W www.dutchblue.com

Dutchblue.com Deutschland

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Dutchblue.com United Kingdom

Park House
111 Uxbridge Road
Ealing, London, W5 5TL

T 00 800 2611 2611 (free) or 02 033 977 242
F 02 033 977 524
E info@dutchblue.com
W www.dutchblue.com

dutchblue.com

DROPPIES

Gebruikershandleiding – Valmat

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User Manual – Fall mat

REF 4100
IFU 4100 – 13M2

dutchblue.com

DROPPIES fall mat

The Dutchblue DROPPIES bedside fall mat offers simple and effective protection in situations where nursing home resident risk a fall from the bed. DROPPIES bedside fall mats are made from high-impact foam with an easy-to-clean vinyl cover. The non-slip material ensures that the floor mat remains firmly in place. The sloping edges make it easy to manoeuvre a wheelchair over the DROPPIES.

Instructions for use

- Remove the packaging. Do not use scissors, knives or other sharp implements when opening the packaging.
- Remove the floor mat from the carry bag and unfold.
- Place the mat on the floor next to the nursing bed. The sloping edges will be visible.
- To limit the drop we recommend lowering the nursing bed as close to the floor as possible.

Instructions for cleaning

The floor mat can be cleaned with a soft cloth and mild detergents or clinical cleaning products. Never use disinfection agents, phenol or solvents.

PLEASE NOTE:

- Do not allow to come into contact with open fire.
- Avoid contact with sharp objects.

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