Dutchblue.com Nederland

Markt 30a
2611 GV Delft
Nederland

T 00 800 2611 2611 (gratis) of 015 3030 025
F 084 883 66 06
E info@dutchblue.com
W www.dutchblue.com

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Square de Meeûs 37
1000 Brussel
België

T 00 800 2611 2611 (gratis) of (0032) 026 699 663
F (0032) 02 403 05 55
E info@dutchblue.com
W www.dutchblue.com

Dutchblue.com Deutschland

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Dutchblue.com United Kingdom

Park House
111 Uxbridge Road
Ealing, London, W5 5TL

T 00 800 2611 2611 (free) or 02 033 977 242
F 02 033 977 524
E info@dutchblue.com
W www.dutchblue.com

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User Manual – Anti-decubitus seat pad

REF 3000
IFU 3000 – 13M2

SITTIES Anti-decubitus seat pad

Sit comfortably on one of our SITTIES anti-decubitus seat pads. The seating pad moulds itself to the contours of the body and thus distributes pressure as evenly as possible. Ideal for the prevention and cure of decubitus complaints.

Instructions for use

Inflating the seat pad

- Remove the packaging. Do not use scissors, knives or other sharp implements when opening the packaging.
- Open the valve by turning it anti-clockwise three times.
- The seat pad will inflate automatically. This may take a few minutes.
- Once the seat pad has been inflated properly, turn the valve clockwise to shut.
- Place the seat pad on the chair or wheelchair.

Deflating the seat pad

- To deflate the seat pad, open the valve completely.
- Press the air out of the seat pad, working in the direction of the valve.
- You will hear the air coming out of the seat pad.
- Close the valve carefully once the seat pad has been fully deflated.

Instructions for cleaning

The seat pad can be cleaned with a soft cloth and mild detergent or clinical cleaning products. Never use disinfection agents, phenol or solvents.

PLEASE NOTE:

- The typical odour of visco-elastic foam will disappear in time.
- Do not allow to come into contact with open fire.
- Avoid contact with sharp objects.

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